


# BLU 200 + C 200 + S 500 + ECP 12

Gentle processing of sensitive parenterals. Seamless, integrated, and highly flexible.



## Secondary packaging of sensitive parenteral products

The BLU 200 blister machine and C 200 cartoner: two perfectly coordinated components for entry into the highly-flexible secondary packaging of parenteral products. With a capacity of up to 200 blisters and 200 cartons per minute, the line also packs syringes and vials in combination blisters made of easily recyclable materials such as PP-Mono or PET. Packing of the cartons into cases and stacking on pallets is automated with the ECP 12 end-of-line packer, the innovative combination of case packer and palletizing module.

Available with medical device tracking, serialization optional directly in the cartoner or with the S 500 as well as aggregation information about the entire packaging line.



REGISTER NOW

### Learn more about the packaging line for parenterals

Register now for the **virtual Uhlmann Plaza** and find out about our packaging line for parenterals.

### Output data at a glance:

<b>BLU 200</b>	max. 200 blisters/min.
<b>C 200</b>	max. 120/200* cartons/min.
<b>S 500</b>	max. 300/500* cartons/min.
<b>ECP 12</b>	max. 12 cases/min. max. 6 pallet cycles/min.

\* optional speed upgrade

